



### Experiment 16 Bragg Reflection


Name: .....

#### Aim

To measure the distance between the lattice planes in a crystal by using X-ray diffraction.

#### Set-up

The set-up consists of a rotatable crystal, placed between a fixed X-ray tube and a Geiger-Müller tube. The GM tube is connected to a rate meter.



The GM tube is mounted on a rotatable measuring arm. The rotations of the crystal and the GM tube are coupled: to rotate the crystal over an angle  $\alpha$ , the GM tube has to be rotated over an angle of  $2 \cdot \alpha$ . As a consequence, the GM-tube always measures the intensity of the X rays reflected by the lattice planes of the crystal (see the diagram above).

The X-ray tube produces X rays with two different wavelengths: 0.138 nm and 0.154 nm. When varying the angle  $\alpha$ , the X rays with a wavelength  $\lambda$  of 0.154 nm give the highest maximum of reflected radiation.

Read the introduction on page 19 of the booklet *ISP Experiments* about X-ray diffraction.

#### Measurements

- Switch on the X-ray tube: set the *time interval* button to 50 minutes, switch on the *power* button, and press the X-ray button. The device is ready for measuring when the red indicator lights up.
- Move the measuring arm with the GM tube along the scale of the device between roughly  $25^\circ$  and  $50^\circ$ . By observing the rate meter, find the angle where the count rate clearly increases quite fast for a slight variation of the angle. The angle found in this way is roughly the first order maximum for the 0.154 nm wavelength radiation.  
Angle (roughly): ..... $^\circ$
- Position the measuring arm at an angle  $1^\circ$  less than the angle found above, and measure the intensity  $I$  (in pulses per 10 s) of the reflected radiation two times. Shift the measuring arm over an angle of  $0.5^\circ$ , and repeat the intensity measurements. Continue to do this until the measuring arm is at an angle of  $1^\circ$  past the angle found in measurement task 2. Record your measurements in the table below. For each of the angles, calculate the average intensity  $I_{avr}$  of the reflected radiation, and record your results in the table below.

observed angle GM tube ( $2 \cdot \alpha$ )	angle between beam and crystal ( $\alpha$ )	$I$ (pulses/10s)		$I_{avr}$ (pulses/10s)

**Note:** In working out the assignments below, you will need the equipment set-up twice for checking your answers in Assignments 3 and 4.

### Assignments

- 1 The angle when measuring the highest intensity is the first order maximum. Calculate the distance  $d$  between the lattice planes (or: the lattice constant) of the crystal. In doing this, use the formula  $n \cdot \lambda = 2 \cdot d \cdot \sin \alpha$  (see the booklet *ISP Experiments*, page 19).

.....  
.....

Lattice constant of the crystal:  $d = \dots\dots\dots$  nm

- 2 The table below shows the lattice constants of some substances. Identify the substance of the crystal in this experiment.

substance	LiF	NaCl	KCl	RbCl
lattice constant $d$ (nm)	0.201	0.281	0.314	0.327

The crystal used is .....

- 3 The X-ray device emits radiation with a wavelength  $\lambda$  of 0.154 nm. This radiation will also give a second order maximum. Calculate the angle  $\alpha$  of this second order maximum. In doing this, use the lattice constant  $d$  found in Assignment 1 and the formula  $n \cdot \lambda = 2 \cdot d \cdot \sin \alpha$ .

.....  
.....

Second order maximum for  $\lambda = 0.154$  nm:  $\alpha = \dots\dots^\circ$

Check your answer by measuring this angle with the equipment set-up.

- 4 The X-ray device also emits radiation with a wavelength  $\lambda$  of 0.138 nm. Also this radiation will give a first order maximum. Calculate the angle  $\alpha$  of this first order maximum.

.....  
.....

First order maximum for  $\lambda = 0.138$  nm:  $\alpha = \dots\dots^\circ$

Check your answer by measuring this angle with the equipment set-up.