

Experiment 18 Elastic Modulus of Rubber


Name:

Aim

To measure the elastic modulus of rubber.

Set-up

The set-up consists of a strip of rubber, clamped to a force meter and a spanner (for stretching the strip of rubber). A source of strontium-90 (⁹⁰Sr) and a Geiger-Müller tube are mounted on opposite sides of the strip of rubber. The GM tube is connected to a pulse counter.



Read the introduction on page 22 of the booklet *ISP Experiments* about the way to measure the elastic modulus.

Measurements



With the equipment set-up the elastic modulus E of rubber can be determined by measuring the force F on and the elongation $\Delta\ell$ of the strip of rubber, in combination with the cross-section area A and the original length ℓ_0 of the strip. Determining the cross-section area A requires measuring the width w and thickness d of the strip of rubber. This thickness d is measured through the absorption of β -radiation in the strip of rubber. The measured intensity of the transmitted radiation and the calibration graph on the other side of this worksheet allow the reading of the thickness d . All other quantities can be measured directly by means of a force meter and a ruler.

- 1 Measure the length ℓ_0 of the strip of rubber (without stretching) between the fixed end points.
 $\ell_0 = \dots\dots\dots$ cm
- 2 Set the force meter to 20 N with help of the spanner. Measure the length ℓ and the width w of the strip of rubber. Record your measurements in the table below. Measure the intensity I (in pulses per 50 s) of the transmitted radiation. Do this by measuring the number of pulses in 10 s five times without resetting the counter. Record your measurement in the table below.
- 3 Repeat this for the other values of the force F as indicated in the table.

| F (N) | ℓ (cm) | w (cm) | $\Delta\ell$ (cm) | I (pulses/50s) | I_{avr} (pulses/10s) | d (cm) | A (cm ²) | $A \cdot \Delta\ell$ (cm ³) |
|------------|----------------|-------------|----------------------|---------------------|---------------------------|-------------|---------------------------|--|
| 20 | | | | | | | | |
| 15 | | | | | | | | |
| 10 | | | | | | | | |
| 5 | | | | | | | | |

Assignments

- 1 For each of the values of the force F on the strip of rubber, calculate the elongation $\Delta\ell$ ($\Delta\ell = \ell - \ell_0$), and record the results in the table above.
- 2 For each of the values of the force F on the strip of rubber, calculate the average intensity I_{avr} of the transmitted radiation in pulses per 10s. Then use the calibration graph at the other side of this worksheet to determine the corresponding thickness d of the strip of rubber. And finally, calculate the corresponding values of the cross-section area A of the strip of rubber and the product of A and $\Delta\ell$. Record all results in the table above.
- 3 Plot the force F as a function of the product $A \cdot \Delta\ell$ in the graph on the other side of this worksheet (right).



- 4 As you can see from the graph, the quantities F en $A \cdot \Delta \ell$ are proportional. Or, in other words: the value of $F/(A \cdot \Delta \ell)$ is a constant. From the graph, determine the value of this constant.

$$\frac{F}{A \cdot \Delta \ell} = \dots\dots\dots \text{ N/cm}^3$$

- 5 Calculate the elastic modulus E of rubber by using the answer of Assignment 4 and the formula that gives the relationship between the force F and the elongation $\Delta \ell$ of the strip of rubber (see the booklet *ISP Experiments*, page 22):

$$F = \frac{E \cdot A \cdot \Delta \ell}{l_0}$$

Hint: First write this formula in the form $E = \dots$

$$E_{\text{rubber}} = \dots\dots\dots \text{ N/cm}^2$$